

Ethnicity: The Region IV Chlamydia project requires patient ethnicity data using the following codes:

- 1 White, Non Hispanic
 - 2 Black, Non Hispanic
 - 4 American Indian
 - 5 Asian
 - 6 Hispanic White
 - 7 Hispanic Black
-

Use the Nucleic Acid Amplification Test (NAAT) transport provided, after collection store at 2° - 30°C

Submissions should occur within 24 hours of collection.

Urine specimens must be transferred to the urine transport tube within 24 hours of collection.

Patient's name and date of collection must be included on specimen container.

Specimen collection kits are supplied by this laboratory.

Only use specimen collection kits within the stamped expiration date.